

Recombinant Mouse Nesfatin-1

CATALOG #:	4913-10
AMOUNT:	10 µg
LOT #:	_____
MOL. WT.	8 kDa
SOURCE:	Produced in <i>E. coli</i> (a.a. 25-106 of mouse nesfatin-1)
PURITY:	>95% by SDS-PAGE
ENDOTOXIN:	<0.1 EU/µg purified protein (LAL test)
FORM:	Lyophilized from PBS

RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute with sterile water to a concentration of 0.1-1 mg/ml. Further dilutions should be made with medium containing 0.2% BSA.

STORAGE CONDITIONS:

The lyophilized protein should be stored desiccated below 0°C. Stable for 6 months. Reconstituted Nesfatin-1 should be stored in working aliquots at -20°C. Avoid freeze/thaw cycles.

DESCRIPTION:

Nesfatin-1, a polypeptide produced by cleavage of the NEFA/nucleobindin-2 (NUCB2), is a secreted protein, expressed in the appetite-control hypothalamic nuclei. Nesfatin is involved in the regulation of the feeding behavior. Intracerebroventricular injection of nesfatin-1 reduces body weight in rats affecting the animal food intake.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies

FOR RESEARCH USE ONLY! Not to be used on humans.