

Recombinant Murine VEGF-120

CATALOG #:	4964-10	10 µg
	4964-100	100 µg
	4964-1000	1 mg
SOURCE:	<i>E. coli</i>	
PURITY:	> 95% by SDS-PAGE	
ENDOTOXIN LEVELS:	Endotoxin level is ≤ 1 EU/µg of rm-VEGF-120.	
MOLECULAR WT.:	Dimer, 14.2/28.4 kDa (with 121/242 amino acids)	

FORMULATION: Lyophilized from a sterile (0.2 micron) filtered aqueous solution containing 0.1% Trifluoroacetic Acid (TFA)

RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in sterile water to a concentration of 0.1mg/ml. The solution can then be diluted into other aqueous buffers and stored at 4°C for 1 week or -20°C for future use.

STORAGE CONDITIONS:

The lyophilized recombinant murine VEGF-120 is very stable at -20°C. Reconstituted rm-VEGF-120 should be stored at working aliquots at -20°C. It is recommended to add a carrier protein (0.1 % HSA or BSA) for long term storage.

DESCRIPTION:

Vascular Endothelial Growth factor-A (VEGF-A) was originally isolated from tumor cells and is produced by a wide variety of cell types. In addition to stimulating vascular growth and permeability, VEGF-A may play a role in stimulating vasodilatation via nitric oxide-dependent pathways. Mouse VEGF-A has several variants, one being VEGF-120 (also known as VEGF-2). Recombinant Mouse VEGF-120 is a homodimer with a molecular weight of 28.4 kDa.

BIOLOGICAL ACTIVITY:

The activity as determined by dose-dependent proliferation of human umbilical vein endothelial cells (HUVEC) was typically found to be 1-5 ng/ml ($\geq 2.0 \times 10^5$ units/mg)

AA SEQUENCE:

MAPTTTEGEQKSHEVIKFMVDVYQRSYCRPIETLVDIFQEYYPDEIEYIFKPCVPLMRCAGCCND
EALCECVPTSESINITMQIMRIKPHQSQHIGEMSFLLQHSRCECRPKKDRTKPEKCDKPRR

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies

FOR RESEARCH USE ONLY! Not to be used in humans.