

Anti-human BCA-1 Polyclonal Antibody

CATALOG NO: 5002-30T 30 µg (Trial size)
5002-100 100 µg

HOST: Rabbit

IMMUNOGEN: *E. coli* expressed human BCA-1 (Internal ID# BV-A56)


FORMULATION:
0.5 mg/ml affinity purified rabbit anti-human BCA polyclonal antibody in phosphate-buffered saline (PBS) containing 30% glycerol, 0.5% BSA, and 0.01% thimerosal.

SPECIES REACTIVITY: human

STORAGE CONDITIONS:
Store at -20°C. For long term storage, aliquot and freeze at -70°C. Avoid repeated freeze/thaw cycles.

BACKGROUND:
Human BCA-1 (B Cell-Attracting chemokine 1), is also known as BLC or CXCL13. It is a recently identified new CXC chemokines. Human BCA-1 is a highly effective attractant for human blood B lymphocytes. The antibody recognizes a 10.1 kDa of recombinant human BCA-1 protein (Cat. #:4001-100)

APPLICATION:
The antibody can be used in Western blotting (0.5-4 µg/ml), However, the optimal concentrations should be determined individually.



Cat. #:5002-100
Western blot analysis of human BCA-1 using recombinant human BCA-1.
Lane 1: 50 ng
Lane 2: 250 ng
Lane 3: 1 µg

FOR RESEARCH USE ONLY! Not to be used on humans.

RELATED PRODUCTS:

- Apoptosis Detection Kits & Reagents
 - Annexin V Kits & Bulk Reagents
 - Caspase Assay Kits & Reagents
 - Mitochondrial Apoptosis Kits & Reagents
 - Nuclear Apoptosis Kits & Reagents
 - Apoptosis Inducers and Set
 - Apoptosis siRNA Vectors
- Cell Fractionation System
 - Mitochondria/Cytosol Fractionation Kit
 - Nuclear/Cytosol Fractionation Kit
 - Membrane Protein Extraction Kit
 - Cytosol/Particulate Rapid Separation Kit
 - Mammalian Cell Extraction Kit
 - FractionPREP Fractionation System
- Cell Proliferation & Senescence
 - Quick Cell Proliferation Assay Kit
 - Senescence Detection Kit
 - High Throughput Apoptosis/Cell Viability Assay Kits
 - LDH-Cytotoxicity Assay Kit
 - Bioluminescence Cytotoxicity Assay Kit
 - Live/Dead Cell Staining Kit
- Cell Damage & Repair
 - HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
 - HAT Colorimetric Assay Kit & Reagents
 - DNA Damage Quantification Kit
 - Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits
- Signal Transduction
 - cAMP & cGMP Assay Kits
 - Akt & JNK Activity Assay Kits
 - Beta-Secretase Activity Assay Kit
- Adipocyte & Lipid Transfer
 - Recombinant Adiponectin, Survivin, & Leptin
 - CETP Activity Assay & Drug Discovery Kits
 - PLTP Activity Assay & Drug Discovery Kits
 - Total Cholesterol Quantification Kit
- Molecular Biology & Reporter Assays
 - siRNA Vectors
 - Cloning Insert Quick Screening Kit
 - Mitochondrial & Genomic DNA Isolation Kits
 - 5 Minutes DNA Ligation Kit
 - 20 Minutes Gel Staining/Destaining Kit
 - β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit
- Growth Factors and Cytokines
- Monoclonal and Polyclonal Antibodies