

Anti-Human IL-17 Polyclonal Antibody

CATALOG #:	5176-100
AMOUNT:	100 µg
LOT #:	_____
IMMUNOGEN:	<i>E. coli</i> expressed human IL-17
HOST:	Rabbit

PREPARATION:

Produced from sera of rabbits pre-immunized with highly pure (>98%) recombinant human IL-17. The antibody was purified by affinity chromatography employing immobilized human IL-17 matrix.

FORMULATION:

100 µg (0.5 mg/ml) affinity purified rabbit anti-human IL-17 polyclonal antibody in phosphate buffered saline (PBS), pH 7.2, containing 50% glycerol, 10 mM EDTA, and 0.02% thymersal.

STORAGE CONDITIONS:

Store at -20°C. For long term storage, aliquot and freeze at -70°C. Avoid repeated freeze/defrost cycles.

BACKGROUND DESCRIPTION:

IL-17 was identified from a CD4+ T cells DNA library. IL-17 can be induced from primary peripheral blood CD4+ T cells upon stimulation. The cytokine exhibits a high degree of amino acid identity with HVS13, an open reading frame from a T lymphotropic Herpesvirus saimiri and with murine CTLA8. Supernatant from cells transfected with IL-17 induced IL-16 and IL-18 production and enhanced the surface expression of the intracellular adhesion molecule-1 (ICAM-1) in human fibroblasts. IL-17 is a disulfide-linked homodimer of two 31.0 kDa subunits, each containing 136 amino acid residues.

APPLICATION AND USAGE:

The antibody can be used for Western blot analysis (0.5-1 µg/ml) and ELISA (2-4 µg/ml). However, the optimal conditions should be determined individually.

FOR RESEARCH USE ONLY! Not to be used on humans.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay Kit
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Polyclonal and Monoclonal Antibodies