


## IL-21 Polyclonal Antibody

<b>CATALOG NO:</b>	5182-30T      30 µg (Trial size) 5182-100      100 µg
<b>CONCENTRATION:</b>	0.5 mg/ml
<b>IMMUNOGEN:</b>	<i>E. coli</i> expressed recombinant human IL-21 (Cat. No. 4182-10)
<b>INTERNAL ID:</b>	BV- C69
<b>HOST:</b>	Rabbit
<b>SPECIES REACTIVITY:</b>	Human
<b>PURIFICATION:</b>	Affinity purified rabbit IgG
<b>FORMULATION:</b>	Supplied in phosphate buffered saline (PBS), pH 7.2, containing 30% glycerol, 0.5% BSA, 0.01% thimerosal
<b>STORAGE CONDITIONS:</b>	Store at -20°C. For long term storage, aliquot and freeze at -70°C. Avoid repeated freeze/thaw cycles.
<b>DESCRIPTION:</b>	IL-21 is a homodimeric, non-glycosylated polypeptide containing 133 amino acid residues and has a molecular weight of 15.5 kDa. IL-21 is a cytokine produced by activated T-cells and can elicit the proliferation of T and B cells and natural killer cells (NK). From recent research, IL-21 has the potential to serve as a target for immunotherapy treatment for treating diseases.
<b>APPLICATION:</b>	Western blot: 0.5-4 µg/ml

**Note: This information is only intended as a guide. The optimal dilutions must be determined by the user.**



**Western blot analysis of IL-21 using recombinant human IL-21**

**IL-21**

### RELATED PRODUCTS:

#### Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

#### Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

#### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

#### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

#### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

#### Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

#### Growth Factors and Cytokines

#### Monoclonal and Polyclonal Antibodies

**FOR RESEARCH USE ONLY! Not to be used on humans.**