

IFN-alpha 2b, Human CellExp™, human recombinant

CATALOG #:	6459-10	10 µg
	6459-50	50 µg
ALTERNATE NAMES:	Leukocyte interferon, B cell interferon, Type I interferon, IFNA2, IFN-α 2a.	
SOURCE:	Human 293 Cell Expressed	
PURITY:	> 95% by SDS - PAGE	
MOL. WEIGHT:	16 kDa, monomer, glycosylated	
ENDOTOXIN LEVEL:	< 1.0 EU per 1 µg of protein	
FORMULATION:	Lyophilized from a PBS solution.	
RECONSTITUTION:	Reconstitute in sterile PBS containing 0.1% endotoxin-free recombinant human serum albumin.	

STORAGE CONDITIONS:

Aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.

ADVANTAGES:


- Animal-derived product free
- High Activity
- Authentic Glycosylation

DESCRIPTION:

Interferons (IFNs) are proteins made and released by host cells in response to the presence of pathogens. They belong to the large class of glycoproteins known as cytokines. IFN-alpha is produced by macrophages and has antiviral activities. IFNs also have other functions: they activate immune cells, such as natural killer cells and macrophages; they increase recognition of infection or tumor cells by up-regulating antigen presentation to T lymphocytes; and they increase the ability of uninfected host cells to resist new infection by virus.

BIOLOGICAL ACTIVITY:

ED₅₀ is typically 0.02 to 0.08 ng/mL. The specific activity was determined by the dose-dependent cytotoxicity of the human TF-1 cell line (human erythroleukemic indicator cell line).



Human Cell^{exp} Human Recombinant IFN-α 2b

RELATED PRODUCTS:

- Human Cell^{exp} Human Recombinant IFN-alpha 2a (Cat # 6458-10, -50)
- IFN-alpha 2a, human recombinant (Cat # 4594-20, -100, -1000)
- IFN-alpha 2b, human recombinant (Cat # 4595-20, -100, -1000)
- IFN-alpha 1, human recombinant (Cat # 4596-10, -100, -1000)
- IFN-beta 1b, human recombinant (Cat # 4860-10, -50, -1000)
- IFN-gamma, human recombinant (Cat # 4116-20, -100, -1000)
- IFN-gamma, murine recombinant (Cat # 4117-20, -100, -1000)
- IFN-gamma, rat recombinant (Cat # 4118-50, -100, -1000)
- IFN-gamma Antibody (Cat # 5116-200)
- IFN-gamma Antibody (Cat # 5117-100)
- IFN-alpha Antibody (Cat # 5596-100)
- IFN-alpha Antibody (Cat # 5595R-100)
- IFN-alpha 2b Antibody (Cat # 5595-100)

FOR RESEARCH USE ONLY! Not to be used in humans.