

Recombinant Protein L

CATALOG #:	6530-1	1 mg
	6530-5	5 mg
	6530-10	10 mg
	6530-100	100 mg
	6530-1000	1 g
	6530- 5000	5 g
LOT #:	_____	
SOURCE:	<i>E. coli</i>	
SEQUENCE:	Amino acid 98-462 of the Peptostreptococcus Magnus. Protein L Contains 5 kappa-binding domain with 6x His-tag on N-terminus. Gene Bank Accession Number AAA25612.	
PURITY:	>98% by SDS-PAGE. Endotoxin level is <0.1 EU/mg of Protein L.	
MOL. WEIGHT:	43.2 kDa with 5 IgG binding sites	
FORM:	Lyophilized with no additives	
RECONSTITUTION:	Reconstitution in water to a concentration of 5 mg/ml gives a clear solution.	

STORAGE CONDITIONS:

Store at -20°C. After reconstitution, aliquot and store at -20°C. Avoid repeated freezing and thawing. Stable, as supplied, for at least 1 year.

APPLICATIONS:

Protein L has the unique ability to bind through kappa light chain interactions without interfering with the antibody's antigen-binding site. This gives Protein L the ability to bind a wider range of Ig classes and subclasses than other antibody-binding proteins. Protein L can be used to detect, quantify and purify antibodies and antibody/antigen complexes. BioVision's Recombinant Protein L contains only IgG binding domains. The albumin-binding domain as well as cell wall and cell membrane binding domains have been removed to ensure the maximum specific IgG binding capacity. The 6-His-tag on N-terminus can be used for affinity purification or protein L detection using anti-His-tag antibodies.

Protein L binds to all classes of Ig (IgG, IgM, IgA, IgE, and IgD). Protein L will also bind Single Chain Variable Fragments (ScFv) and Fab Fragments. Protein L binds kappa I, III, and IV in human and kappa I on mouse.

SPECIFICITY:

Under optimal conditions, 1 mg protein L will bind approximately 5 mg human IgG. Optimal binding of Protein L to antibodies occurs at pH 7.5 over a pH range of 2.5 to 3.0.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay Kit
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Quality Antibodies for Apoptosis and Signal Transduction Molecules

FOR RESEARCH USE ONLY! Not to be used on humans.