


TIP60 Antibody

ALTERNATE NAMES:	KAT5 Antibody: lysine acetyltransferase 5, K-acetyltransferase 5, lysine acetyltransferase 5, HIV-1 Tat interactive protein 60kDa, HTATIP, HTATIP1, TIP, cPLA2 interacting protein, cPLA2, ESA1, PLIP
CATALOG #:	6634-50
AMOUNT:	50 µg
HOST:	Rabbit
IMMUNOGEN:	TIP60 antibody was raised against a 15 amino acid peptide near the carboxy terminus of human TIP60.
PURIFICATION:	Affinity chromatography purified via peptide column
MOLECULAR WEIGHT:	Predicted: 51, 54, 56, 60 kDa. Observed: 53 kDa
FORM:	Liquid
FORMULATION:	Supplied as 1 mg/ml in PBS containing 0.02% sodium azide
SPECIES REACTIVITY:	Human, Mouse, Rat
STORAGE CONDITIONS:	TIP60 antibody can be stored at 4°C for three months and at -20°C, stable for up to one year.


DESCRIPTION: TIP60 belongs to the MYST family of histone acetyl transferases (HATs) and was originally isolated as an HIV-1 TAT interactive protein. HATs play important roles in regulating chromatin remodeling, transcription and other nuclear processes by acetylating histone and non-histone proteins. TIP60 has been shown to play a role in DNA repair and apoptosis and is thought to play an important role in signal transduction and autophagy. In case of HIV-1 infection, interaction with the viral Tat protein leads to TIP60 polyubiquitination and targets it to degradation. Also, TIP60 contains a chromodomain. TIP60 is a catalytic subunit of the NuA4 histone acetyltransferase complex. The NuA4 HAT complex plays a role in transcriptional activation of select genes mainly by acetylation of nucleosomal histone H4 and H2A, which influence nucleosome-DNA interaction and promotes interaction of the modified histones with other proteins that could regulate transcription positively. This complex also involves in the activation of transcriptional programs associated with oncogene and proto-oncogene mediated growth induction, tumor suppressor mediated growth arrest and replicative senescence, DNA repair and apoptosis.

APPLICATION: Western blot: 1 - 2 µg/mL, IHC: 2.5 µg/mL, IF: 20 µg/ml.


Note: This information is only intended as a guide. The optimal dilutions must be determined by the user.



Western blot analysis of TIP60 in human brain tissue lysate with TIP60 antibody at 1 µg/mL.



Immunohistochemistry of TIP60 in rat brain tissue with TIP60 antibody at 2.5 µg/ml.



Immunofluorescence of TIP60 in rat brain tissue with TIP60 antibody at 20 µg/ml.

RELATED PRODUCTS:

- TIP60 Antibody (Cat. No. 6126-50)
- HAT-1 Antibody (Cat. No. 3689-100)
- HAT-1 Blocking Peptide (Cat. No. 3689BP-50)
- HAT-2 Antibody (Cat. No. 3692-50)
- HAT-2 Blocking Peptide (Cat. No. 3692BP-50)
- HAT-3 Antibody (Cat. No. 3707-50)
- HAT-3 Blocking Peptide (Cat. No. 3707BP-50)

FOR RESEARCH USE ONLY! Not to be used on humans.