


SOD1 Antibody

CATALOG NO:	6641-30T 30 µg (Trial size) 6641-100 100 µg
ALTERNATE NAMES:	Superoxide Dismutase 1, SODC
HOST:	Rabbit
ISOTYPE:	IgG
IMMUNOGEN:	Recombinant Human SOD1
INTERNAL ID:	BV-M61
PURIFICATION:	Affinity purified rabbit IgA
FORM:	Liquid
FORMULATION:	0.5 mg/ml of antibody in PBS pH 7.2, containing 30% glycerol, 0.5% BSA, 5 mM EDTA and 0.03% Proclin.
SPECIES REACTIVITY:	Human, Mouse and Rat.
STORAGE CONDITIONS:	Store for 1 year at -20°C from date of shipment. Avoid repeated freeze/thaw cycles.

DESCRIPTION: Superoxide Dismutase (SOD) is an oxidoreductase that catalyzes the reaction between superoxide anions and hydrogen to yield molecular oxygen and hydrogen peroxide. The enzyme protects the cell against dangerous levels of superoxide. It belongs to the Cu-Zn superoxide dismutase family. It binds copper and zinc ions. The encoded isozyme is a soluble cytoplasmic protein, acting as a homodimer to convert naturally-occurring but harmful superoxide radicals to molecular oxygen and hydrogen peroxide. The other isozyme is a mitochondrial protein. Mutations in this gene have been implicated as causes of familial amyotrophic lateral sclerosis. Rare transcript variants have been reported for this gene. SOD1 destroys radicals which are normally produced within the cells and which are toxic to biological systems. Defects in SOD1 are the cause of amyotrophic lateral sclerosis type 1 (ALS1). ALS1 is a familial form of amyotrophic lateral sclerosis, a neurodegenerative disorder affecting upper and lower motor neurons and resulting in fatal paralysis. Sensory abnormalities are absent. Death usually occurs within 2 to 5 years. The etiology of amyotrophic lateral sclerosis is likely to be multifactorial, involving both genetic and environmental factors. The disease is inherited in 5-10% of cases leading to familial forms.

APPLICATION: Western blot: 1:200

Note: This information is only intended as a guide. The optimal dilutions must be determined by the user.



Western blot with SOD1 antibody

Lane 1: 3T3 cell lysate.
Lane 2: Rat kidney lysate.

RELATED PRODUCTS:

- Superoxide Dismutase, human recombinant (**Cat. No. 4802-20, -100, -1000**)
- Human Recombinant SOD2 (**Cat. No. 6160-100**)
- Bacterial Recombinant SOD-A (**Cat. No. 6361-100**)
- SOD-1 Antibody (**Cat. No. 3458-100**)
- SOD-1 Blocking Peptide (**Cat. No. 3458BP-50**)
- Superoxide Dismutase 1 Antibody (72B1) (**Cat. No. 6168-100**)
- Superoxide Dismutase 2 Antibody (2A1) (**Cat. No. 6169-100**)
- Superoxide Dismutase 3 Antibody (1H12) (**Cat. No. 6170-100**)
- Superoxide Dismutase 4 Antibody (3A1) (**Cat. No. 6171-100**)
- Superoxide Dismutase (SOD) Activity Assay Kit (**Cat. No. K335-100**)

FOR RESEARCH USE ONLY! Not to be used on humans.