

Goat Anti-Rat IgG (H&L) FITC

CATALOG #:	6907-250
AMOUNT:	250 µl
HOST/ISOTYPE:	Goat
IMMUNOGEN:	Rat IgG (H&L)
FORM:	Liquid
FORMULATION:	In 0.01 M Sodium phosphate, 0.25 M NaCl, 50% Glycerol, pH 7.6 and 0.05% Sodium Azide.
PURIFICATION:	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
SPECIES REACTIVITY:	Rat.
STORAGE CONDITIONS:	Upon delivery aliquot and store at -20°C for six months. Avoid freeze/thaw cycles.
DESCRIPTION:	Goat Polyclonal Secondary Antibody to Rat IgG (H&L) FITC.
SPECIFICITY:	By immunoelectrophoresis and ELISA this antibody reacts specifically with rat IgG and with light chains common to other rat immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins.
APPLICATION:	IF: 1:50 – 1:200, FACS: 1:50 – 1:200.

Note: This information is only intended as a guide. The optimal dilutions must be determined by the user.

RELATED PRODUCTS:

- Anti-Rabbit IgG, HRP-Linked Antibody (Cat # 6401-05)
- Anti-Mouse IgG, HRP-Linked Antibody (Cat # 6402-05)
- Anti-Goat IgG, HRP-Linked Antibody (Cat # 6403-05)
- Goat Anti-Mouse IgG (H&L) Cy3 (Cat # 6900-250)
- Goat Anti-Mouse IgG (H&L) Texas Red (Cat # 6901-250)
- Goat Anti-Mouse IgG (H&L) AMCA (Cat # 6902-250)
- Goat Anti-Rabbit IgG (H&L) Cy3 (Cat # 6903-250)
- Goat Anti-Rabbit IgG (H&L) Texas Red (Cat # 6904-250)
- Goat Anti-Rabbit IgG (H&L) AMCA (Cat # 6905-250)
- Goat Anti-Rat IgG (H&L) Cy3 (Cat # 6906-250)
- Goat Anti-Rat IgG (H&L) HRP (Cat # 6908-250)
- Goat Anti-Rat IgG (H&L) TRITC (Cat # 6909-250)
- Goat Anti-Rat IgG (H&L) Biotin (Cat # 6910-250)
- Rabbit Anti-Goat IgG (H&L) Cy3 (Cat # 6911-250)

FOR RESEARCH USE ONLY! Not to be used on humans.