

Human CellExp™ CD304 / NRP1 / Neurophilin-1, human recombinant

CATALOG #: 7375-10 10 µg

ALTERNATE NAMES: NRP1, Neuropilin-1, NRP, VEGF165R, CD304

SOURCE: HEK 293 cells (Phe 22 – Lys 644)

PURITY: ≥ 92% by SDS-PAGE gel

MOL. WEIGHT: This protein is fused with polyhistidine tag at the C-terminus, has a calculated MW of 70.7 kDa. The predicted N-terminus is Phe 22. DTT-reduced Protein migrates as 85-90 kDa due to glycosylation.

ENDOTOXIN LEVEL: <1 EU/µg by LAL method

FORM: Lyophilized


FORMULATION: Lyophilized from 0.22 µm filtered solution in PBS, pH7.4. Generally 5-8% Mannitol or trehalose is added as a protectant before lyophilization.

STORAGE CONDITIONS: Store at -20°C. After reconstitution, aliquot and store at -20°C and use within 3 months. Avoid repeated freezing and thawing cycles.

RECONSTITUTION: Centrifuge the vial prior to opening. Reconstitute in sterile PBS, pH 7.4 to a concentration of 50 µg/ml. Do not vortex. This solution can be stored at 2-8°C for up to 1 month. For extended storage, it is recommended to store at -20°C.

DESCRIPTION: Neuropilin-1 (NRP1) is also known as vascular endothelial cell growth factor 165 receptor (VEGF165R), CD antigen CD304, which belongs to the neuropilin family. The membrane-bound isoform 1 is a receptor involved in the development of the cardiovascular system, in angiogenesis, in the formation of certain neuronal circuits and in organogenesis outside the nervous system. It mediates the chemorepulsant activity of semaphorins. It binds to semaphorin 3A, The PLGF-2 isoform of PGF, The VEGF-165 isoform of VEGF and VEGF-B. Coexpression with KDR results in increased VEGF-165 binding to KDR as well as increased chemotaxis. It may regulate VEGF-induced angiogenesis. The soluble isoform 2 binds VEGF-165 and appears to inhibit its binding to cells.

BIOLOGICAL ACTIVITY: Measured by its binding ability in a functional ELISA. Immobilized rh NRP1 / Neuropilin-1 / CD304 can bind rhVEGF165 with an apparent K_d < 1 nM.



Human recombinant CD304/NRP1/Neurophilin-1

RELATED PRODUCTS:

- Human CellExp™ CD223, human recombinant (Cat. No. 7278-10, -50)
- Human CellExp™ CD71, human recombinant (Cat. No. 7279-10, -50)
- Human CellExp™ CD273, human recombinant (Cat. No. 7369-10, -50)
- Human CellExp™ CD33, human recombinant (Cat. No. 7370-10, -50)
- Human CellExp™ CD36, human recombinant (Cat. No. 7371-10, -50)
- Human CellExp™ CD87, human recombinant (Cat. No. 7372-20, -100)
- Human CellExp™ CD360, human recombinant (Cat. No. 7373-20, -100)
- Human CellExp™ CD244, human recombinant (Cat. No. 7374-10, -50)

FOR RESEARCH USE ONLY! Not to be used in humans.