

# Human CellExp™ VEGF120, mouse recombinant

**CATALOG #:** 7424-10 10 µg

**ALTERNATE NAMES:** VEGF120, VEGFA, VPF

**SOURCE:** HEK 293 cells (Ala 27 - Arg 146)

**PURITY:** ≥ 95% by SDS-PAGE gel

**MOL. WEIGHT:** Mouse VEGF120, contains no "tag", has a calculated MW of 17.2 kDa. The predicted N-terminus is Ala 27. DTT-reduced Protein migrates as 18-22 kDa due to glycosylation.

**ENDOTOXIN LEVEL:** <1 EU/µg by LAL method

**FORM:** Lyophilized


**FORMULATION:** Lyophilized from 0.22 µm filtered solution in PBS. Generally 5-8% Mannitol or trehalose is added as a protectant before lyophilization.

**STORAGE CONDITIONS:** Store at -20°C. After reconstitution, aliquot and store at -20°C and use within 3 months. Avoid repeated freezing and thawing cycles.

**RECONSTITUTION:** Centrifuge the vial prior to opening. Reconstitute in sterile PBS, pH 7.4 to a concentration of 50 µg/ml. Do not vortex. This solution can be stored at 2-8°C for up to 1 month. For extended storage, it is recommended to store at -20°C.

**DESCRIPTION:** Vascular endothelial growth factor A (VEGFA) is also known as Vascular permeability factor (VPF). VEGFA belongs to the PDGF/VEGF growth factor family. VEGFA is a glycosylated mitogen that specifically acts on endothelial cells and has various effects, including mediating increased vascular permeability, inducing angiogenesis, vasculogenesis and endothelial cell growth, promoting cell migration, and inhibiting apoptosis. Alternatively spliced transcript variants, encoding either freely secreted or cell-associated isoforms, have been characterized. VEGFA is produced by a group of three major isoforms as a result of alternative splicing and if any three isoforms are produced (VEGFA120, VEGFA164, and VEGFA188) then this will not result in vessel defects and death of the full VEGFA knockout in mice.

**BIOLOGICAL ACTIVITY:** Measured in a cell proliferation assay using human umbilical vein endothelial cells (HUVEC). The ED<sub>50</sub> for this effect is typically 5 ng/ml.



**RELATED PRODUCTS:**

- Human CellExp™ Human Recombinant VEGF-C (Cat # 7231-10)
- Human CellExp™ Human Recombinant VEGF 165 (Cat # 6485-10, -50)
- VEGF121, human recombinant (Cat. No. 4963-10, -50, -1000)
- VEGF165, human recombinant (Cat. No. 4363-10, -50, -1000)
- VEGF165, murine recombinant (Cat. No. 4364-10, -50, -1000)
- VEGF165, rat recombinant (Cat. No. 4365-10, -50, -1000)
- VEGF120, murine recombinant (Cat. No. 4964-10, -100, -1000)
- VEGF-B, human recombinant (Cat. No. 4642-10, -20, -1000)
- VEGF-C, human recombinant (Cat. No. 4633-10, -50, -1000)
- VEGF-C, murine recombinant (Cat. No. 4634-10, -50, -1000)
- VEGF-C, rat recombinant (Cat. No. 4635-10, -50, -1000)
- VEGF-D, human recombinant (Cat. No. 4343-10, -50, -1000)

**FOR RESEARCH USE ONLY! Not to be used in humans.**