

Human CellExp™ IL7Ra/CD127, human recombinant

CATALOG #: 7492-10 10 µg
7492-50 50 µg

ALTERNATE NAMES: IL7Ra, CD127

SOURCE: HEK 293 cells (Glu 21 – Gly 236)

PURITY: ≥ 95% by SDS-PAGE gel

MOL. WEIGHT: This protein is fused with Fc fragment of human IgG1 at the C-terminus, has a calculated MW of 51 kDa. The predicted N-terminus is Glu 21. DTT-reduced Protein migrates as 66-90 kDa due to glycosylation.

ENDOTOXIN LEVEL: <1 EU/µg by LAL method

FORM: Lyophilized

FORMULATION: Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally Mannitol or Trehalose is added as protectants before lyophilization.


STORAGE CONDITIONS: Store at -20°C. After reconstitution, aliquot and store at -20°C and use within 3 months. Avoid repeated freezing and thawing cycles.

RECONSTITUTION: Centrifuge the vial prior to opening. Reconstitute in 50 mM Tris, 100 mM glycine, pH 7.0. Do not vortex. This solution can be stored at 2-8°C for up to 1 month. For extended storage, it is recommended to store at -20°C.

DESCRIPTION: Interleukin-7 receptor subunit alpha (IL7Ra) also known as CD antigen CD127, which belongs to the type I cytokine receptor family and type 4 subfamily. IL7Ra /CD127 contains 1 fibronectin type-III domain. IL7Ra /CD127 is expressed on various cell types, including naive and memory T cells and many others. IL7Ra /CD127 is receptor for interleukin-7 and also acts as a receptor for thymic stromal lymphopoietin (TSLP). The IL7 receptor is a heterodimer of IL7R and IL2RG. The TSLP receptor is a heterodimer of CRLF2 and IL7R.

BIOLOGICAL ACTIVITY: Measured by its binding ability in a functional ELISA.

Immobilized rhIL7Ra Fc Chimera at 0.25 µg/ml (100 µl/well) on a Rabbit anti-human IgG Fc antibody, the concentration of rhIL-7 that produces 50% of the optimal binding response is found to be approximately 20 - 100 ng/mL.



RELATED PRODUCTS:

- Human CellExp™ CD223, human recombinant (Cat. No. 7278-10, -50)
- Human CellExp™ CD71, human recombinant (Cat. No. 7279-10, -50)
- Human CellExp™ CD273, human recombinant (Cat. No. 7369-10, -50)
- Human CellExp™ CD33, human recombinant (Cat. No. 7370-10, -50)
- Human CellExp™ CD36, human recombinant (Cat. No. 7371-10, -50)
- Human CellExp™ CD87, human recombinant (Cat. No. 7372-20, -100)
- Human CellExp™ CD360, human recombinant (Cat. No. 7373-20, -100)
- Human CellExp™ CD244, human recombinant (Cat. No. 7374-10, -50)
- Human CellExp™ CD304, human recombinant (Cat. No. 7375-10)
- Human CellExp™ CD319, human recombinant (Cat. No. 7376-10, -50)
- Human CellExp™ CD306, human recombinant (Cat. No. 7377-10, -50)
- Human CellExp™ CD84, human recombinant (Cat. No. 7378-10, -50)

FOR RESEARCH USE ONLY! Not to be used in humans.