


## Histone H4 (1-103 aa), Human recombinant

CATALOG #:	7671-100	100 µg
	7671-250	250 µg
ALTERNATE NAMES:	Histone H4, H4.	
SOURCE:	<i>E. coli</i>	
PURITY:	≥95% by SDS - PAGE	
MOL. WEIGHT:	11.5 kDa (1-103 aa)	
FORMULATION:	Supplied in water	
STORAGE CONDITIONS:	Store at -80°C. Avoid repeated freeze and thaw cycles.	

### DESCRIPTION:

A nucleosome is the basic repeating unit of chromatin in which 146 base pairs of DNA wrap twice around a histone octamer consisting of two copies of each of the core histones, H2A, H2B, H3, and H4. The combination of two H2A/H2B dimers and one H3/H4 tetramer create the nucleosome core. Histone H3 undergoes many modifications which include acetylation, methylation, and phosphorylation that are important for regulation of gene transcription.



### RELATED PRODUCTS:

- Histone Protein (**Cat. No. 1136-10**)
- Histone H2A (1-129 aa), Xenopus recombinant (**Cat. No. 7668-100, -250**)
- Histone H2B (1-123 aa), Xenopus recombinant (**Cat. No. 7669-100, -250**)
- Histone H3 (1-136 aa), Human recombinant (**Cat. No. 7670-100, -250**)
- Histone H2A Antibody (**Cat. No. 3621-100**)
- Histone H2AX Antibody (**Cat. No. 3761-100**)
- Histone H2B Antibody (**Cat. No. 3622-100**)
- Histone H3 Antibody (**Cat. No. 3623-100**)
- Histone H3 Peptide-Biotin Conjugate (**Cat. No. 2233-01, -1**)
- Histone H4 Antibody (**Cat. No. 3624-100**)
- Histone Methyltransferase (SUV39H1) Antibody (**Cat. No. 3942-100**)
- Acetyl-Histone H2A Antibody (**Cat. No. 3653-100**)
- Acetyl-Histone H2B Antibody (**Cat. No. 3654-100**)
- Acetyl-Histone H4 Antibody (**Cat. No. 3656-100**)

**FOR RESEARCH USE ONLY! Not to be used in humans.**