

Recombinant Human MMP-1 Catalytic Domain

CATALOG #:	7781-10	10 µg
	7781-50	50 µg
	7781-1000	1 mg
LOT #:	_____	
SYNONYMS:	Matrix Metalloproteinase-1, Interstitial Collagenase, Fibroblast Collagenase	
SOURCE:	<i>E. coli</i>	
MOL. WT.:	19.9 kDa	
PURITY:	>95% by SDS-PAGE and HPLC Endotoxin level is <0.1 ng per µg of MMP-1.	
FORM:	Lyophilized powder	
RECONSTITUTION:	Centrifuge the vial prior to opening. Reconstitute in sterile PBS to a concentration of 0.1-1.0 mg/ml. This solution can then be diluted into other aqueous buffers and stored at 4° C for 1 week or -20° C for future use.	
STORAGE CONDITIONS:	The lyophilized protein is best-stored desiccated below 0° C. Reconstituted MMP1 should be stored in working aliquots at -20° C. For long-term storage, it is recommended to add carrier protein (0.1% HAS or BSA).	
DESCRIPTION:	Active MMP-1 catalytic domain from human cDNA, expressed in <i>E.coli</i> and purified using proprietary technologies. The active MMP-1 protein is very useful in studying enzyme kinetics, cleave target substrates, and screen for inhibitors, etc.	
BIOLOGICAL ACTIVITY:	The protein is biologically active as determined by MMP-1 activity assay.	

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies

FOR RESEARCH USE ONLY! Not to be used on humans.