

ToxOut™ Endotoxin-free Collection Tube, 2.0 ml

rev 11/16

Cat. No.
7937-50 ToxOut™ Endotoxin-free Collection Tube, 2.0 ml, 50 Tubes

I. Description:

Endotoxin is the lipopolysaccharide (LPS) complex located in the outer membrane of gram-negative bacteria. During experimental procedures, large amount of endotoxins are shed and can easily contaminate labware, buffers and downstream products. Presence of endotoxin in any of these has a detrimental effect on culturing of cells and can lead to generation of irreproducible or negative experimental data. Therefore, it is critical to keep endotoxin away from samples and products.

ToxOut™ Endotoxin-free Collection Tube, 2.0 ml are screw cap microcentrifuge tubes that are recommended for samples with special requirement of endotoxin level (< 0.05 EU/ml). BioVision's Endotoxin-free tubes are strong enough for centrifugation up to 20,000 x *g* and suitable for use in vapor-phase liquid nitrogen storage, -80 °C freezers, boiling, and autoclaving.

II. Storage and Handling:

Store the tubes at room temperature.

Handle endotoxin free tube with extra caution to prevent endotoxin contamination from dust, solution or dirty labware.



III. Related Products:

- ToxOut™ Rapid Endotoxin Removal Kit (Cat. No. K2501-5)
- ToxOut™ Rapid Endotoxin Removal Agarose (Cat. No. 7941-5)
- ToxOut™ High Capacity Endotoxin Removal Kit (Cat. No. K2502-5)
- ToxOut™ High Capacity Endotoxin Removal Agarose (Cat. No. 7942-5)
- ToxOut™ Endotoxin Rapid Endotoxin Removal Regeneration Buffer (Cat No. 7939-50)
- ToxOut™ Endotoxin Removal Equilibration Buffer (Cat No. 7940-50)
- ToxOut™ Endotoxin-free Water (Cat. No. 7938-50)
- ToxOut™ Endotoxin-free PBS (Cat. No. 7943-50)
- ToxOut™ Endotoxin-free Glass Vial (Cat. No. 7944-2, -10)

FOR RESEARCH USE ONLY! Not to be used on humans.

155 S. Milpitas Blvd., Milpitas, CA 95035 USA | T: (408)493-1800 F: (408)493-1801 | www.biovision.com | tech@biovision.com