



Anti-MYCBP Antibody

CATALOG NO:	A1199-100
ALTERNATIVE NAMES:	AMY1, C-Myc-binding protein, Associate of Myc 1, AMY-1
AMOUNT:	100 µl
IMMUNOGEN:	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human MYCBP
HOST/ISOTYPE:	Rabbit IgG
CLONALITY:	Polyclonal
SPECIFICITY:	Recognizes endogenous levels of MYCBP protein
SPECIES REACTIVITY:	Human, Mouse and Rat
PURIFICATION:	The antibody was purified by affinity chromatography
FORM:	Liquid
FORMULATION:	Supplied in 0.42% Potassium phosphate; 0.87% Sodium chloride; pH 7.3; 30% glycerol and 0.01% sodium azide
STORAGE CONDITIONS:	Shipped at 4°C. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
DESCRIPTION:	May control the transcriptional activity of MYC. Stimulates the activation of E box-dependent transcription by MYC
APPLICATION:	WB; 1:500 – 1:2000, IHC; 1:50 – 1:200

Note: This information is only intended as a guide. The optimal dilutions must be determined by the user.



Immunohistochemical analysis of MYCBP staining in human heart formalin fixed paraffin embedded tissue section



Western blot analysis of MYCBP expression in A549 (A); rat lung (B) whole cell lysates

RELATED PRODUCTS:

- MBD1 polyclonal antibody (Cat. No. 6828-25)
- MCAF Antibody (Cat. No. 3759-100)
- Anti-ZNF7 Antibody (Cat. No. A1163-100)
- NIF1 Antibody (Cat. No. 3737-100)

FOR RESEARCH USE ONLY! Not to be used on humans.