

UltraPolymer Goat Anti-Rabbit/Mouse IgG (H&L) HRP Cocktail

CATALOG NO:	A1280-5 5 ml A1280-15 15 ml
HOST/ISOTYPE:	Goat
IMMUNOGEN:	Purified Rabbit IgG, whole molecule; Purified Mouse IgG, whole molecule
FORM:	Liquid
FORMULATION:	In PBS with 1% BSA and 0.1% proclin 150
PURIFICATION:	Affinity purified using solid phase Sheep IgG
SPECIES REACTIVITY:	Rabbit, Mouse
STORAGE CONDITIONS:	Store at 2-8°C.
DESCRIPTION:	Universal UltraPolymer Cocktail, Goat anti-Rabbit/Mouse IgG (H&L) conjugated to HRP, affinity purified

SPECIFICITY: Based on IEP, the Goat anti-Rabbit IgG immunopolymer antibody reacts with heavy (γ) chains on rabbit IgG, light chains on all rabbit immunoglobulins; the Goat anti-Mouse IgG immunopolymer antibody reacts with heavy (γ) chains on mouse IgG, light chains on all mouse immunoglobulins;
Cross-Reactivity: Based on IEP, no cross-reactivity is observed from Goat-anti-Rabbit IgG to non-immunoglobulin rabbit serum proteins, serum proteins from human, mouse or rat; no cross-reactivity is observed from Goat-anti-Mouse IgG to non-immunoglobulin mouse serum proteins, serum proteins from bovine, horse, human, pig or rabbit

APPLICATION AND USAGE: IHC: 125- 150 μ l

Note: This information is only intended as a guide. The optimal dilutions must be determined by the user.

RELATED PRODUCTS:

- UltraPolymer Goat Anti-Rabbit IgG (H&L) HRP (**Cat. No. A1274-5**)
- UltraPolymer Goat Anti-Mouse IgG (H&L) HRP (**Cat. No. A1275-5**)
- UltraPolymer Donkey Anti-Goat IgG (H&L) HRP (**Cat. No. A1276-5**)
- UltraPolymer Goat Anti-Rat IgG (H&L) HRP (**Cat. No. A1277-5**)
- UltraPolymer Goat Anti-Human IgG (H&L) HRP (**Cat. No. A1278-5**)
- UltraPolymer Donkey Anti-Sheep IgG (H&L) HRP (**Cat. No. A1279-5**)
- Goat Anti-Rabbit IgG (H&L) Texas Red (**Cat. No. 6904-250**)
- Goat Anti-Rabbit IgG (H&L) AMCA (**Cat. No. 6905-250**)
- Goat Anti-Rat IgG (H&L) Cy3 (**Cat. No. 6906-250**)
- Goat Anti-Rat IgG (H&L) FITC (**Cat. No. 6907-250**)
- Goat Anti-Rat IgG (H&L) HRP (**Cat. No. 6908-250**)
- Goat Anti-Rat IgG (H&L) TRITC (**Cat. No. 6909-250**)
- Goat Anti-Rat IgG (H&L) Biotin (**Cat. No. 6910-250**)
- Rabbit Anti-Goat IgG (H&L) Cy3 (**Cat. No. 6911-250**)

FOR RESEARCH USE ONLY! Not to be used on humans.