


Anti-HES5 Antibody

CATALOG NO:	A1637-100 100 µl
ALTERNATIVE NAMES:	bHLHb38; Class B basic helix-loop-helix protein 38; Hairy and enhancer of split 5; Hairy/Enhancer of split, Drosophila, homolog of, 5; HES 5; hes family bHLH transcription factor 5; Hes5; HES5 protein; HES5_HUMAN; MGI:104876; Transcription factor HES-5
CONCENTRATION:	1 mg/ml
IMMUNOGEN:	A synthesized peptide
HOST/ISOTYPE:	Rabbit IgG
CLONALITY:	Polyclonal
MOL WEIGHT:	19 kDa
SPECIES REACTIVITY:	Human, Mouse, Rat
PURIFICATION:	Affinity purification
FORM:	Liquid
FORMULATION:	Supplied in phosphate buffered saline, pH 7.4, 150 mM NaCl, 0.02% sodium azide and 50% glycerol
STORAGE CONDITIONS:	For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
DESCRIPTION:	Transcriptional repressor of genes that require a bHLH protein for their transcription.
APPLICATION:	WB 1:500-1:2000 IELISA (peptide) 1:20000-1:40000

Note: This information is only intended as a guide. The optimal dilutions must be determined by the user.



Western blot analysis of 293 whole cell lysates using HES5 antibody

RELATED PRODUCTS:

- GFAP Antibody (Cat. No. 3206)
- Vimentin Antibody (Cat. No. 3634)
- Fibronectin Antibody (Cat. No. 3630)
- FGF-2 Antibody (Cat. No. 5039)
- Anti-CD133 Antibody (Cat. No. A1622)
- Anti-CHRD1 Antibody (Cat. No. A1623)
- Anti-HES3 Antibody (Cat. No. A1636)

FOR RESEARCH USE ONLY! Not to be used on humans.