

Anti-Sall4 Antibody

CATALOG NO.: A1692-100

BACKGROUND DESCRIPTION: Transcription factor with a key role in the maintenance and self-renewal of embryonic and hematopoietic stem cells.

ALTERNATE NAMES: DRRS; HSAL4; Sal like 4 (Drosophila); Sal like 4; Sal like Protein 4; Sal-like protein 4; Sall4; SALL4_HUMAN; Spalt like transcription factor 4; Tex20; Zinc finger protein 797; Zinc finger protein SALL4; ZNF797.

AMOUNT: 100 µl.

HOST/ISOTYPE: Rabbit / IgG.

IMMUNOGEN: A synthesized peptide derived from human Sall4.

MOLECULAR WEIGHT: 140 kDa.

PURIFICATION: Affinity purified.


FORM: Liquid.

FORMULATION: In PBS pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol

SPECIES REACTIVITY: Rat.
Mouse.
Human.

STORAGE CONDITIONS: Store at -20°C; For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

APPLICATIONS AND USAGE: Western Blotting.
ELISA.



Western blot analysis of Sall4 using HuvEc whole lysates.

RELATED PRODUCTS:

Anti-KLF4 Antibody (A1684).
Anti-REX 1 Antibody (A1691).
Anti-LRH1 Antibody (A1688).
Anti-MELK Antibody (A1686).

FOR RESEARCH USE ONLY! Not to be used on humans.