

Anti-SP1 Antibody

CATALOG NO.: A1999-100 (100 µl)

BACKGROUND DESCRIPTION: The protein encoded by this gene is a zinc finger transcription factor that binds to GC-rich motifs of many promoters. The encoded protein is involved in many cellular processes, including cell differentiation, cell growth, apoptosis, immune responses, response to DNA damage, and chromatin remodeling. Post-translational modifications such as phosphorylation, acetylation, glycosylation, and proteolytic processing significantly affect the activity of this protein, which can be an activator or a repressor. Three transcript variants encoding different isoforms have been found for this gene.

ALTERNATE NAMES: TSFP1; Transcription factor Sp1

ANTIBODY TYPE: Polyclonal

HOST/ISOTYPE: Rabbit / IgG

IMMUNOGEN: KLH-conjugated synthetic peptide targeting a sequence within the center region of human SP1

MOLECULAR WEIGHT: 90 kDa

PURIFICATION: Affinity purified

FORM: Liquid


FORMULATION: In 0.42% Potassium phosphate; 0.87% NaCl; pH 7.3; 30% glycerol; and 0.01% sodium azide

SPECIES REACTIVITY: Human, Mouse, Rat, Monkey, Chicken, Bovine


STORAGE CONDITIONS: Store at -20°C. Avoid freeze / thaw cycles

APPLICATIONS AND USAGE: WB 1:500 - 1:1000, IF 1:100 - 1:500

Note: This information is only intended as a guide. The optimal dilutions must be determined by the user



Western blot analysis of SP1 expression in K562 (A); A431 (B) whole cell lysates.



Immunofluorescent analysis of SP1 staining in K562 cells. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 5-10 minutes and blocked with 3% BSA-PBS for 30 minutes at RT. Cells were probed with the primary antibody in 3% BSA-PBS and incubated overnight at 4 C in a humidified chamber. Cells were washed with PBST and incubated with a DyLight 594-conjugated secondary antibody (red) in PBS at RT in the dark. DAPI was used to stain the cell nuclei (blue).

RELATED PRODUCTS:

Anti-p53 Antibody (A1805)

Anti-STAT3 Antibody (A1670)

Phospho SP1 (Thr453) Antibody (A1979)

Smad4 Antibody (3464)

FOR RESEARCH USE ONLY! Not to be used on humans.