


Anti-NRF2 Antibody

CATALOG NO.: **A2019-100 (100 µl)**


BACKGROUND DESCRIPTION: This gene encodes a transcription factor which is a member of a small family of basic leucine zipper (bZIP) proteins. The encoded transcription factor regulates genes which contain antioxidant response elements (ARE) in their promoters; many of these genes encode proteins involved in response to injury and inflammation which includes the production of free radicals. Multiple transcript variants encoding different isoforms have been characterized for this gene.

ALTERNATE NAMES:	NRF2; Nuclear factor erythroid 2-related factor 2; NF-E2-related factor 2; NFE2-related factor 2; HEBP1; Nuclear factor erythroid derived 2 like 2
ANTIBODY TYPE:	Polyclonal
HOST/ISOTYPE:	Rabbit / IgG
IMMUNOGEN:	KLH-conjugated synthetic peptide targeting a sequence within the C-term region of human NRF2
MOLECULAR WEIGHT:	68 kDa
PURIFICATION:	Affinity purified
FORM:	Liquid
FORMULATION:	In 0.42% Potassium phosphate; 0.87% NaCl; pH 7.3; 30% glycerol; and 0.01% sodium azide
SPECIES REACTIVITY:	Human, Mouse, Rat, Porcine
STORAGE CONDITIONS:	Store at -20°C. Avoid freeze / thaw cycles
APPLICATIONS AND USAGE:	WB 1:500 - 1:1000, IHC 1:100 - 1:200, IF 1:100 - 1:500


Note: This information is only intended as a guide. The optimal dilutions must be determined by the user



Western blot analysis of NRF2 expression in THP1 (A); HL60 (B); K562 (C); Jurkat (D); mouse kidney (E); rat kidney (F) whole cell lysates.



Immunohistochemical analysis of NRF2 staining in human lung cancer formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0), then incubated with the antibody at RT and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with hematoxylin and mounted with DPX.



Immunofluorescent analysis of NRF2 staining in HL60 cells. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 5-10 minutes and blocked with 3% BSA-PBS for 30 minutes at RT. Cells were probed with the primary antibody in 3% BSA-PBS and incubated overnight at 4 C in a humidified chamber. Cells were washed with PBST and incubated with a DyLight 594-conjugated secondary antibody (red) in PBS at RT in the dark. DAPI was used to stain the cell nuclei (blue).

RELATED PRODUCTS:

GSK-3b Antibody (3494)
AKT1 Antibody (NT) (6745)
HDAC2 Antibody (3602)
Anti-Catalase Antibody (A1311)

FOR RESEARCH USE ONLY! Not to be used on humans.