

Anti-IL-6R α , Human IgG1 Antibody (Sr2)

rev 11/21

CATALOG NO.: A2133-100 (100 μ g)

BACKGROUND DESCRIPTION: The research-grade biosimilar is a fully human anti-IL-6R monoclonal IgG1 antibody that can bind to membrane as well as soluble forms of interleukin 6 receptor (IL-6R).**ALTERNATE NAMES:** IL6Q, gp80, CD126, IL6RA, IL6RQ, IL-6RA, IL-6R-1, HIES5, IL-6R, IL-1Ra, IL6QTL, IL-6R-Alpha, IL6R, IL-6R 1**ANTIBODY TYPE:** Monoclonal**HOST/ISOTYPE:** Recombinant / Human IgG1**SOURCE:** CHO cells**IMMUNOGEN:** Human IL-6 receptor alpha**FORM:** Liquid**FORMULATION:** In PBS buffer, pH 7.2**STORAGE CONDITIONS:** Undiluted antibody solution should be stored at 4 °C in the dark and avoid freezing the solution**APPLICATIONS AND USAGE:** In vitro and In vivo studies**This information is only intended as a guide. The optimal dilutions must be determined by the user****RELATED PRODUCTS:**

Anti-IL-6 receptor, Human IgG1 Antibody (A1447)

Anti-PD-1, Humanized Antibody (A1306)

Anti-IL-6 Antibody (A1398)

Anti-PD-L1, Human IgG1 Antibody (A2109)

Anti-IL-2R alpha (CD25), Human IgG1 Antibody (A1098)

FOR RESEARCH USE ONLY! Not to be used on humans.