

Anti-SARS-CoV-2 RBD Antibody (Clone# CB6)

CATALOG NO.: A2286-50 (50 µg)
A2286-100 (100 µg)

BACKGROUND DESCRIPTION: Coronaviruses (CoV) are a large group of enveloped positive-sense RNA viruses. They belong to subfamily Coronavirinae, in the family of Coronaviridae, of the order of Nidovirales. The Coronavirus genome is about 30 kb in length and encodes four structural proteins, namely, spike (S), envelope (E), membrane (M) and nucleocapsid (N), multiple nonstructural proteins, and other accessory proteins. The Spike (S) protein is a 180 kDa type I glycosylated transmembrane protein, that assembles as a trimer on the surface of the virus, hence the coronavirus appears crown-shaped (In Latin, corona means crown). The ectodomain of the Spike protein consists of two domains: the S1 domain whose main function is receptor binding, and the S2 domain whose main function is membrane fusion. S1 binds to the cell surface receptor of the host for attachment during the viral entry, this induces a conformational change in the S2, thus enabling the fusion of host and viral membranes and ultimately entry of the viral genome in the host cell.

ALTERNATE NAMES: Spike protein, COVID19, COVID 19, S protein, SARS-CoV S protein, S glycoprotein, E2, Peplomer protein, Spike protein S1, SARS Coronavirus, SARS-CoV-2, SARS CoV 2, 2019-nCoV, RBD Protein, CB6

ANTIBODY TYPE: Monoclonal

CLONE: CB6

HOST/ISOTYPE: Recombinant / Human IgG

IMMUNOGEN: Recombinant SARS-CoV-2 RBD protein

FORM: Liquid


FORMULATION: In PBS, pH 7.5

SPECIES REACTIVITY: Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), SARS-CoV


STORAGE CONDITIONS: Store at -20 to -80°C. Avoid repeated freeze/thaw cycles

APPLICATIONS AND USAGE: ELISA 1:5,000 - 1:10,000

This information is only intended as a guide. The optimal dilutions must be determined by the user



SDS-PAGE analysis was performed to assess the purity and integrity of Anti-SARS-CoV-2 RBD Antibody (Clone# CB6).



Direct ELISA shows that Anti-SARS-CoV-2 RBD Antibody (Clone# CB6) can bind to immobilized RBD (EC₅₀ = 72.47 ng/ml).

RELATED PRODUCTS:

Anti-SARS-CoV-2 Spike (NTD) Antibody (Clone# 4A8) (A2269)
 Anti-CoV-2 & SARS-CoV S1 Antibody (Clone# CR3022) (A2103)
 Anti-SARS-CoV-2 Spike S1 Antibody (A3000)
 Anti-SARS-CoV-2 Spike S1 Antibody (Clone# 4C6) (A3001)
 Anti-SARS-CoV-2 S1 Antibody (Clone# S309) (A2266)

FOR RESEARCH USE ONLY! Not to be used on humans.