


PRODUCT: GSK2292767

ALTERNATE NAMES: GSK 2292767; GSK-2292767; N-(5-(4-(5-(((2R,6S)-2,6-dimethylmorpholino)methyl)oxazol-2-yl)-1H-indazol-6-yl)-2-methoxy pyridin-3-yl)methanesulfonamide

CATALOG#: B1660-5, -25

AMOUNT: 5 mg, 25 mg

STRUCTURE:



MOLECULAR FORMULA: C₂₄H₂₈N₆O₅S

MOLECULAR WEIGHT: 512.58

CAS NUMBER: 1254036-66-2

APPEARANCE: Solid powder

SOLUBILITY: DMSO

PURITY: >98%

STORAGE: Dry, dark and at 4°C for short term (days to weeks) or at -20°C for long term. Protect from air and light.

DESCRIPTION: GSK2292767 is a potent and selective PI3K δ inhibitor. GSK2292767 is currently in clinical trials for the treatment of respiratory diseases such as asthma and COPD.

HANDLING: Do not take internally. Wear gloves and mask when handling the product! Avoid contact by all modes of exposure.

REFERENCES: 1. Down KD et al. Optimization of Novel Indazoles as Highly Potent and Selective Inhibitors of Phosphoinositide-3-Kinase Delta for the Treatment of Respiratory Disease. J Med Chem. 2015 Aug 24.

RELATED PRODUCTS:

- APY-0201 (Cat. No. 9564-5,25)
- AS-604850 (Cat. No. 2339-1, 5)
- BAG956 (Cat.No. 2456-5, 25)
- BAY 80-6946 (Cat. No. B1267-5,25)
- BEZ235 (NVP-BEZ235) (Cat. No. 1626-5, 25)
- CAL-101 (Cat. No. 2613-5, 25)
- GDC-0941 bismesylate (Cat. No. 1623-1, 5)
- Deforolimus (Cat. No. 1587-5, 25)
- DiscoveryPak™ PI 3-Kinase Inhibitor Panel (Cat. No. K856-5)
- GDC-0941 (Cat. No. 1623-1,5)
- GSK-263771 (Cat. No. 2793-5, 25)
- IC87114 (Cat. No. 1661-1)
- IPI-145 (Cat. No. 2445-5, 25)
- LY 294002 (Cat. No. 1667-5, 25)
- NVP-BKM120 (Cat. No. 2619-5, 25)
- mTOR Inhibitor, Ku-0063794 (Cat. No. 1779-1,5)
- OSU-03012, hydrochloride (Cat. No. 2295-5, 25)
- PI-103 (Cat. No. 1728-1,5)
- PIK-294 (Cat. No. 2614-5, 25)
- PI3-Ky Inhibitor, AS-605240 (Cat. No. 1780-1,5)
- PathwayReady™ PI3-K/Akt/mTOR Signaling Inhibitor panel (Cat. No. K857-11)
- PP242 (Cat. No. 1658-1)
- PX-866 (Cat. No. 1965-1, 5)
- Rapamycin (Cat. No. 1568-5, 50)
- SAR260301 (Cat. No. B1263-5,25)
- Temsirolimus (Cat. No. 1600-5, 25)
- TGX-115 (Cat. No. 1660-1)
- TGX-221 (Cat. No. 1781-1,5)
- Wortmannin (Cat. No. 1670-1)
- YM201636 (Cat. No. 2045-1, 5)
- XL147 (Cat. No. 9401-5, 25)
- TGX-115 (Cat. No. 1660-1)
- TGX-221 (Cat. No. 1781-1, 5)
- Wortmannin (Cat. No. 1670-1)

USAGE: *FOR RESEARCH USE ONLY! Not to be used in humans*