


PRODUCT: **EZSolution™ BMS-299897,
Sterile-Filtered**

ALTERNATE NAME (R)-4-(2-(1-(4-chloro-N-(2,5-difluorophenyl)phenylsulfonamido)ethyl)-5-fluorophenyl)butanoic acid

CATALOG #: B2243-5

AMOUNT: 5 mg

STRUCTURE:



MOLECULAR FORMULA: C₂₄H₂₁ClF₃NO₄S

MOLECULAR WEIGHT: 511.94

CAS NUMBER: 290315-45-6

APPEARANCE: Liquid

FORMULATION: A 10 mM sterile-filtered solution in DMSO

PURITY: ≥98%

STORAGE: Store at -20°C. Protect from air and light

DESCRIPTION: A 10 mM (5 mg in 977 µl), sterile-filtered solution of the γ-secretase inhibitor BMS-299897 (Cat. No. 9606) in anhydrous DMSO.

RELATED PRODUCTS:

APP Antibody (Cat. No. 6001-100)
Beta-Secretase, human recombinant (Cat. No. 7604-20, 100, 1000)
Beta-Secretase Activity Assay Kit (Cat. No. K360-100)
Beta Secretase Inhibitor I (Cat. No. 7501-1)
Beta Secretase Inhibitor II (Cat. No. 7502-1)
BMS-288897 (Cat. No. 9606-5,25)
EZSolution™ BMS-299897 (Cat. No. B2243)
Compound 34 (Cat. No. 1898-200, 1000)
Compound E (Cat. No. 1949-250, 1000)
Compound W (Cat. No. 2208-50, 250)
DAPT (Cat. No. 1855-5)
EZSolution™ DAPT (Cat. No. 2663-5)
DBZ (Cat. No. 2649-1,5)
L-685, 258 (Cat. No. 1713-250)
LY2886721 (Cat. No. 2399-5, 25)
LY-411575 (Cat. No. 2648-1, 5)
RO4929097 (Cat. No. 2011-1, 5)
Semagacestat (Cat. No. 2430-1, 5)
TACE Antibody (Cat. No. 3222-100)

FOR RESEARCH USE ONLY! Not to be used on humans.

HANDLING: Do not take internally. Wear gloves and mask when handling the product! Avoid contact by all modes of exposure.