


EZSolution™ Peficitinib

ALTERNATE NAME: ASP015K
JNJ-54781532
4-(((1R,2r,3S,5s,7s)-5-hydroxyadamantan-2-yl)amino)-1H-pyrrolo[2,3-b]pyridine-5-carboxamide
JNJ54781532

CATALOG #: B2623-1 1 mg

STRUCTURE:



MOLECULAR FORMULA: C₁₈H₂₂N₄O₂

MOLECULAR WEIGHT: 326.4

CAS NUMBER: 944118-01-8

APPEARANCE: Liquid

PURITY: ≥98% by HPLC

DESCRIPTION: A 10 mM (1 mg in 306 µl) of the JAK inhibitor Peficitinib (Cat. No. B2474) in anhydrous DMSO

STORAGE TEMPERATURE: -20°C. Protect from moisture

HANDLING: Do not take internally. Wear gloves and mask when handling the product! Avoid contact by all modes of exposure.

RELATED PRODUCTS:

EZSolution™ PF-03394197 (B2479)
EZSolution™ Baricitinib (B2480)
EZSolution™ Cyt387 (B2484)
Peficitinib (B2474)

DISCLAIMER: *FOR RESEARCH USE ONLY! Not to be used on humans.*