


Toltrazuril

ALTERNATE NAME: 1-methyl-3-[3-methyl-4-[4-[(trifluoromethyl)thio]phenoxy]phenyl]-1,3,5-triazine-2,4,6(1H,3H,5H)-trione
BAY-Vi-9142

CATALOG #: B3404-1000 1 g
B3404-250 250 mg

STRUCTURE:



MOLECULAR FORMULA: C₁₈H₁₄F₃N₃O₄S

MOLECULAR WEIGHT: 425.4

CAS NUMBER: 69004-03-1

APPEARANCE: Solid

PURITY: ≥98%

SOLUBILITY: 25 mg/ml DMSO

DESCRIPTION: Toltrazuril is an antiprotozoal agent that is widely used to prevent coccidia parasite infection in veterinary applications.

STORAGE TEMPERATURE: -20°C. Protect from light

HANDLING: Do not take internally. Wear gloves and mask when handling the product! Avoid contact by all modes of exposure.

RELATED PRODUCTS:

Atovaquone (B2640)
Curvularin (9694)
Sinefungin (B2824)
Furazolidone (B2697)

DISCLAIMER: ***FOR RESEARCH USE ONLY! Not to be used on humans.***