

Screw Caps

(Cat# M1353-100; -250; -500; Store at RT)

I. Introduction:

BioVision's Screw Caps are molded of high-purity, low metal-content polypropylene copolymer resin. The Screw caps are pressure-tested together at 7.5 psig for safe air shipment.

II. Features:

- Low Profile, nonpyrogenic, noncytotoxic
- Vials sold separately (BV# M1352)
- Nonsterile
- Leakproof when used with BioVision's Vials (BV# M1352)
- Color: White
- Closure size 13 mm
- **This item is recommended for BioVision's Vials (BV# M1352) in Cat# K206, K226, K406 and K555**

III. Package Contents:

Cat. No.	Quantity/Pack
M1353-100	100
M1353-250	250
M1353-500	500

IV. Shipment and Storage:

The Screw Caps should be stored at RT.

V. Related Products:

Cat. No.	Product Name
M1352-100	Polypropylene Vials
M1352-250	Polypropylene Vials
M1352-500	Polypropylene Vials

FOR RESEARCH USE ONLY! Not to be used on humans.