

96-Well Black Plates (2)

CATALOG #: M1355-2

AMOUNT: 2 plates

INTRODUCTION:

These are polystyrene microplates (solid black) with standard solid flat bottom, chimney well. The flat bottom allows precise optics and the rounded corners allow pipetting without leaving a residue.


FEATURES:

- Free of detectable DNase, RNase, human DNA
- Non-pyrogenic
- Packaging dimension: 335 x 250 x 165 mm
- The risk of sample carryover and cross-contamination is minimized
- **Suitable for fluorescence measurements**
- The polystyrene surface is more hydrophobic than the surface of a high binding microplate, and therefore tends to be more suitable for non polar proteins and peptides
- Solid bottom
- Without lids
- Autoclavability: No

STORAGE CONDITIONS: Room temperature

APPLICATIONS:

These plates are useful for diagnostic and immunological research applications and ideal for fluorometric assays using a 96-well format and standard plate readers.



RELATED PRODUCTS:

- 10K Spin Column
- 96-Well Plate-2 (transparent)
- 96-Well Plate-2 (white)
- Avidin Covalinked 96-well Plates
- Dounce Tissue Grinder
- Magnetic Separator (Re-usable)
- Protein G Coated 96-well Plates
- Metabolic Assay Kits
- Caspase Assay Kits

FOR RESEARCH USE ONLY! Not to be used on humans.