


Isopropyl Alcohol

ALTERNATE NAME: Isopropanol, IPA, sec-Propyl alcohol, 2-Propanol

CATALOG #: M1526-5 5 ml
M1526-50 50 ml



MOLECULAR FORMULA: (CH₃)₂CHOH

MOLECULAR WEIGHT: 60.10

CAS NUMBER: 67-63-0

APPEARANCE: Clear liquid

PURITY: ≥ 99%

DESCRIPTION: Isopropanol is the simplest of the secondary alcohols. It is flammable and miscible in chloroform, water, ethanol, and ether. It has applications as a disinfectant, a solvent, a sample preservative etc.

STORAGE TEMPERATURE: Store at room temperature

HANDLING: Do not take internally. Wear gloves and mask when handling the product! Avoid contact by all modes of exposure.

RELATED PRODUCTS:

K808 (Deproteinizing Sample Preparation Kit)
K555 (Hydroxyproline Colorimetric Assay Kit)
K823 (Deproteinizing Sample Preparation Kit II)

DISCLAIMER: ***FOR RESEARCH USE ONLY! Not to be used on humans.***