

# CA Syringe Filter, 0.22 $\mu$ m, 25 mm

(Luer-Lok/Luer Slip, Sterile 100 or 200/pack, Store at RT)

<b>CATALOG #:</b>	M4201-100 M4201-200
<b>QUANTITY:</b>	100 per pack 200 per pack
<b>MEMBRANE:</b>	Cellulose Acetate
<b>FILTRATION AREA(CM2):</b>	4.08
<b>WETTABILITY:</b>	Hydrophilic
<b>MAXIMUM OPERATING TEMPERATURE:</b>	50°C
<b>BURST PRESSURE (PSI):</b>	87
<b>BUBBLE POINT (PSI):</b>	0.28
<b>FLOW RATE (ml/min@10psi):</b>	24
<b>HOLDUP VOLUME:</b>	100
<b>HOUSING MATERIAL:</b>	PP
<b>STORAGE:</b>	Room Temperature

**DESCRIPTION:** Cellulose Acetate Syringe Filters are low protein binding syringe filters. Hence these filters are ideal for tissue culture media and sensitive biological samples filtration. Filters also contain Luer-Lok inlet/Luer Slip connector for a secure fit to the syringe. These are individually packed to maintain the sterility of the filters.



Cellulose Acetate Syringe Filter, Sterile

## RELATED PRODUCTS:

- Cellulose Acetate Syringe Filters (Cat# M4200-M4205).
- PES Syringe Filters (Cat# M4302- M4306).
- PVDF Syringe Filters (Cat# M4469- M4472)

**FOR RESEARCH USE ONLY! Not to be used on humans.**