

CA Syringe Filter, 0.45 μ m, 33 mm

(Luer-Lok/Luer Slip, Sterile, 100 or 200/pack, Store at RT)

CATALOG #:	M4205-100 M4205-200
QUANTITY:	100 per pack 200 per pack
MEMBRANE:	Cellulose Acetate
FILTRATION AREA(CM²):	4.9
WETTABILITY:	Hydrophilic
MAXIMUM OPERATING TEMPERATURE:	50°C
BURST PRESSURE (PSI):	87
BUBBLE POINT (PSI):	0.28
FLOW RATE (ml/min@10psi):	70
HOLDUP VOLUME:	100
VOLUME THROUGHPUT(ml):	120
HOUSING MATERIAL:	Polypropylene
STORAGE:	Room Temperature
DESCRIPTION:	Cellulose Acetate Syringe Filters are low protein binding syringe filters. Hence these filters are ideal for tissue culture media and sensitive biological samples filtration. Filters also contain Luer-Lok inlet/Luer Slip connector for a secure fit to the syringe. These are individually packed to maintain the sterility of the filters.



Cellulose Acetate Syringe Filter, Sterile

RELATED PRODUCTS:

- Cellulose Acetate Syringe Filters (Cat# M4200-M4205).
- PES Syringe Filters (Cat# M4302- M4306).
- PVDF Syringe Filters (Cat# M4469- M4472)

FOR RESEARCH USE ONLY! Not to be used on humans.