

CA Syringe Filter, 0.45 μ m, 20 mm, Nonsterile

(Cellulose Acetate filter, Double Luer Lock, 100/pack, Store at RT)

CATALOG #:	M4223-100
QUANTITY:	100 per pack
MEMBRANE:	Cellulose Acetate
HOUSING MATERIAL:	Acrylonitrile butadiene styrene
STORAGE:	Room Temperature

DESCRIPTION: Cellulose Acetate Syringe Filters are low protein binding syringe filters. Hence these filters are ideal for tissue culture media and sensitive biological samples filtration. Filters also contain Double Luer Lock connector for a secure fit to the syringe. These are nonsterile filters that are available in a resealable storage container to prevent contamination during multiple opening and closings.



Cellulose Acetate Syringe Filter, Nonsterile

RELATED PRODUCTS:

- Cellulose Acetate Syringe Filters (Cat# M4216-M4234)
- Glass Fiber Syringe Filters (Cat# M4238-M4248)
- Nylon Syringe Filters (Cat# M4268-M4301)
- PES Syringe Filters (Cat# M4322- M4349)
- Polypropylene Syringe Filters (Cat# M4363- M4384)
- PTFE Syringe Filters (Cat# M4416-M4468)
- PVDF Syringe Filters (Cat# M4483- M4513)

FOR RESEARCH USE ONLY! Not to be used on humans.