

Human CellExp™ Tim-4, (Fc-Tag), Human Recombinant

CATALOG NO:	P1361-10 P1361-50	10 µg 50 µg
ALTERNATE NAMES:	TIM4; TIMD4; T Cell Immunoglobulin and Mucin Domain-containing Protein 4; T Cell Membrane Protein 4	
SOURCE:	HEK 293 cells (aa 25-315)	
PURITY:	≥ 80% by SDS – PAGE	
SEQUENCE:	The extracellular domain of human Tim-4 (aa 25-315) is fused to the N-terminus of the Fc region of mouse IgG2a.	
MOLECULAR WT.:	~80kDa	
ENDOTOXIN LEVEL:	< 5 EU/µg protein by LAL method	
FORM:	Lyophilized	
FORMULATION:	Lyophilized from 0.2 µm-filtered solution in PBS	
STORAGE CONDITIONS:	For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.	
RECONSTITUTION:	Reconstitute with 100 µl sterile water. Add 1X PBS to the desired protein concentration.	
DESCRIPTION:	<p>The TIM (T cell/transmembrane, immunoglobulin and mucin) family plays a critical role in regulating immune responses, including allergy, asthma, transplant tolerance, autoimmunity and the response to viral infections. The unique structure of TIM immunoglobulin variable region domains allows highly specific recognition of phosphatidylserine (PtdSer), exposed on the surface of apoptotic cells. TIM-4 (T cell; immunoglobulin; Mucin-4), also known as SMUCKLER, is a 60 kDa member of the TIM family of immune regulating proteins. TIM-4 is exclusively expressed on antigen-presenting cells, where it mediates phagocytosis of apoptotic cells and plays an important role in maintaining tolerance. TIM-4 binds specifically to TIM-1 which is also the cellular receptor for the hepatitis A virus, and has been implicated in the development of asthma. Among hematopoietic cells, TIM-1 is expressed on activated B and T cells, preferentially in the Th2 subset of CD4+ T cells. The interaction of TIM-4 with TIM-1 induces costimulatory and hyperproliferative signals in T cells.</p>	
BIOLOGICAL ACTIVITY:	Measured by its ability to inhibit anti-CD3-induced proliferation of stimulated human T cells	

RELATED PRODUCTS:

- Human CellExp™ HAVCR1 / KIM1 / TIM1, human recombinant (**Cat. No. 7494**)
- TIM-4, human recombinant (**Cat. No. P1108**)
- Tim-1, Human Recombinant (**Cat. No. P1351**)
- Tim-1, Human Recombinant (**Cat. No. P1352**)
- Human CellExp™ Tim-1, Fc Tag, Human Recombinant (**Cat. No. P1353**)
- Tim-1, Mouse Recombinant (**Cat. No. P1354**)
- Human CellExp™ Tim-2, Fc Tag, Mouse Recombinant (**Cat. No. P1355**)
- Tim-2, Mouse Recombinant (**Cat. No. P1356**)
- Tim-3, Human Recombinant (**Cat. No. P1357**)
- Tim-3, Human Recombinant, (Fc-Tag Mouse) (**Cat. No. P1358**)
- Human CellExp™ Tim-3, Mouse Recombinant (**Cat. No. P1359**)
- Human CellExp™ Tim-3, (Fc-Tag Mouse), Mouse Recombinant (**Cat. No. P1360**)
- Human CellExp™ Tim-4, (Fc-Tag), Mouse Recombinant (**Cat. No. P1362**)
- Tim-4, Mouse Recombinant (**Cat. No. P1363**)

FOR RESEARCH USE ONLY! Not to be used on humans.