

CA 72-4, Human Metastatic Liver Carcinoma

06/19

CATALOG NO: P1436-1 1 KU
P1436-5 5 KU

ALTERNATE NAMES: Cancer Antigen 72-4 protein, CA72-4, CA 72-4, CA72-4 protein

FORM: Liquid

FORMULATION: In 50 mM Phosphate and 150 mM NaCl (pH 7.2 -7.6) containing 0.09% Sodium Azide.

SPECIFIC ACTIVITY: ≥ 10000 U/ml

STORAGE CONDITIONS: Aliquot and store at 4°C. 4.0 °C

DESCRIPTION: Cancer antigen CA 72-4 (CA 72-4) is a biomarker for ovarian cancer. The elevated level of CA 72-4 has been reported in gastrointestinal (GI) cancer patients. Level of Cancer Antigen 72-4 (CA 72-4) has been also been shown to be increased in multiple cancers e.g pancreas, stomach, gallbladder, biliary tract, colon, breast, ovary, cervix, and endometrium. CA 72-4 along with CA 19-9 and CEA assays is used to diagnose and monitor gastrointestinal cancer.

RELATED PRODUCTS:

Leucine-rich Alpha 2 Glycoprotein 1 (LRG1), Human Plasma (P1431)
D-Dimer, Human Plasma (P1434)
Beta 2 Microglobulin, Human Urine (P1432)
Cancer Antigen (CA 15-3), human breast adenocarcinoma (P1435)

FOR RESEARCH USE ONLY! Not to be used on humans.