


## Choose your DNA Polymerase!

DNA Polymerase	BV Cat. No	Extension Speed	Target Length	Proofreading	Specificity as compared to Taq DNA Polymerase (*)	Fidelity as compared to Taq DNA Polymerase (1X)	Special Feature
ExpressTaq™ DNA Polymerase	M1504	1 kb/15 sec	10-15 kb	No	***	1X	Fast and Robust
HiFidelity™ DNA Polymerase	M1505	1kb/20-30 sec	15-20 kb	Yes	****	1300X	Fidelity and Sensitivity
FireTaq™ DNA Polymerase	M1506	1 kb/15 sec	10-15 kb	No	****	1X	High Specificity
Advance™ DNA Polymerase	M1151	1 kb/min	High processivity	Yes	***	50X	For AT/GC rich template

 Best in the category