



# ELISA Kits for Cancer Research

Enzyme-linked Immunosorbent Assay (ELISA) is a powerful analytics tool to quantify proteins and small molecules in serum, plasma, tissue homogenates and other biological fluids. BioVision offers fully validated, accurate, and ready-to-use ELISA Kits for convenient and sensitive quantitation of various cancer and tumor related targets

## Key Features:

- **Accurate:** reproducible results with low intra & inter assay variability
- **High Sensitivity and Specificity**
- **Versatile:** broad range of sample types such as serum, plasma, tissue, cell culture supernatant, urine and many other biological fluids.



**Figure:** Standard curve for human PD-L1 ELISA Kits (Left) and human Interleukin- 10 ELISA kits (right).

Product Name	Cat. No.	Size	Specificity	Detection
3-Nitrotyrosine (3-NT) ELISA Kit	K4158	100 assays	Universal	Competitive
8-hydroxy-2'-deoxyguanosine (8-OHdG) ELISA Kit	K4160	100 assays	Universal	Competitive
AKT1/PKB ELISA Kit	K4166	100 assays	Human	Sandwich
Angiogenin ELISA Kit	K4802	100 assays	Human	Sandwich
Angiopoietin-1 ELISA Kit	K7115	100 assays	Human	Sandwich
Angiopoietin-2 ELISA Kit	K7116	100 assays	Human	Sandwich
Angiostatin ELISA Kit	K4755	100 assays	Human	Sandwich
ANGPTL3 Serum ELISA Kit	K4914	100 assays	Human	Sandwich
ANGPTL3 Serum ELISA Kit	K4915	100 assays	Mouse, Rat	Sandwich
ANGPTL6 Serum ELISA Kit	K4916	100 assays	Human	Sandwich
BDNF ELISA Kit	K4788	100 assays	Human	Sandwich

155 S. Milipitas Blvd, Milipitas, CA 95035

T: 408-493-1800 F: 408-493-1801

Toll Free: 800-891-9699 (US Only)



Product Name	Cat. No.	Size	Specificity	Detection
Beta-NGF ELISA Kit	K4787	100 assays	Human	Sandwich
BLC/BCA-1/CXCL13 ELISA Kit	K4766	100 assays	Human	Sandwich
BMP-4 ELISA Kit	K4745	100 assays	Human	Sandwich
CA19-9 ELISA Kit	K7427	100 assays	Human	Sandwich
Cancer Antigen 125 ELISA Kit	K4803	100 assays	Human	Sandwich
Cancer Antigen 15-3 (CA15-3) ELISA Kit	K4804	100 assays	Human	Sandwich
Carcinoembryonic Antigen (CEA) ELISA Kit	K4805	100 assays	Human	Sandwich
Caspase-9 ELISA Kit	K4165	100 assays	Human	Sandwich
Cathepsin B ELISA Kit	K7408	100 assays	Human	Sandwich
Cathepsin D ELISA Kit	K7409	100 assays	Human	Sandwich
Chorionic Gonadotropin (hCG) ELISA Kit	K7424	100 assays	Human	Sandwich
Cortisol (human/mouse/rat) ELISA Kit	K7430	100 assays	Human, Mouse, Rat	Sandwich
CTLA4 /CD152 ELISA Kit	K4157	100 assays	Human	Sandwich
EGFR/ErbB ELISA Kit	K4168	100 assays	Human	Sandwich
GM-CSF ELISA Kit	K4772	100 assays	Human	Sandwich
Growth Hormone (human) ELISA Kit	K7412	100 assays	Human	Sandwich
HGF ELISA Kit	K4781	100 assays	Human	Sandwich
IL-10 ELISA Kit	K4167	100 assays	Human	Sandwich
IL-1 $\alpha$ ELISA Kit (For Lysates)	K4792	100 assays	Human	Sandwich
IL-1 $\beta$ ELISA Kit	K4794	100 assays	Human	Sandwich
IL-1 $\beta$ ELISA Kit	K4795	100 assays	Mouse	Sandwich
IL-4 ELISA Kit	K4164	100 assays	Human	Sandwich
IL-5 ELISA Kit	K4156	100 assays	Human	Sandwich
IL-7 ELISA Kit	K4154	100 assays	Human	Sandwich
Leptin ELISA Kit	K4777	100 assays	Human	Sandwich
p53 ELISA Kit	K4829	100 assays	Human	Sandwich
PD-1 ELISA Kit	K4153	100 assays	Human	Sandwich
PD-L1 ELISA Kit	K4155	100 assays	Human	Sandwich
Progranulin ELISA Kit	K4738	100 assays	Human	Sandwich
Prostate Specific Antigen ELISA Kit	K7431	100 assays	Human	Sandwich
TGF- $\alpha$ ELISA Kit	K4339	100 assays	Human	Sandwich
TGF- $\beta$ 1 ELISA Kit	K4342	100 assays	Human	Sandwich
TNF- $\alpha$ ELISA Kit	K4779	100 assays	Human	Sandwich
TNF- $\alpha$ ELISA Kit	K1051	100 assays	Mouse	Sandwich
TNF- $\beta$ ELISA Kit	K4780	100 assays	Human	Sandwich
VEGF ELISA Kit	K5363	100 assays	Human	Sandwich

Please visit [www.BioVision.com](http://www.BioVision.com) for a comprehensive overview on these products!