


# PCR, QPCR, and RT PCR Products



BioVision offers a wide range of high performance PCR, QPCR, RT PCR and QRT PCR enzyme products including kits and master mixes, enzymes including DNA Polymerases and Reverse Transcriptases for all your PCR needs. Over the years, our scientists have engineered the highest quality enzymes and most advanced and sophisticated formulations for PCR related applications. All of the BioVision's products have been well validated by researchers around the world, with an increasing number of citations available in scientific publications. BioVision is dedicated to continue contributing to the scientific research community by delivering market leading performing PCR, RTPCR, QPCR products and DNA Polymerases at an unbeatable price.

## Product Highlights:

- Ready to use master mixes
- Highly sensitive, efficient and cost-effective master mixes
- Variety of master mixes to choose from
- Master mixes with dye for direct loading of the PCR product(s) onto an agarose gel
- Wide range of RT PCR enzyme products including kits, master mixes and Reverse Transcriptases
- Numerous DNA Polymerases to choose from
- For routine PCR amplification to PCR applications demanding high efficiency, high yield, high fidelity and robustness
- For amplification of difficult templates

Product Name	Cat. No.	Size
DNA Ladders	9301; M1155-M1156	Varies
DNA Polymerases	9001; 9003; 9222; 9223; 9225; M1146-M1153	Varies
PCR Kits	M1144-M1145	Varies
PCR Master Mixes	M1127-M1143	Varies
QPCR Master Mixes	M1105-M1125	Varies
QPCR Positive Control Kit	M1126	100 Reactions
QRT PCR Products	M1175-M1190	Varies
Reverse Transcriptases	M1173-M1174	Varies
RT PCR Products	M1160-M1172	Varies

We are constantly adding new products to our Molecular Biology Tools category for your research. Please visit our website, <http://www.biovision.com/molecular-biology-tools-979/> for the most up-to-date details on all the products under the Molecular Biology Category.

155 S. Milipitas Blvd, Milipitas, CA 95035

T: 408-493-1800 F: 408-493-1801

Toll Free: 800-891-9699 (US Only)

