


# Secondary Antibodies



Make your **SIGNALS STAND OUT** and get **OUTSTANDING RESULTS**  
with BioVision's Secondary Antibodies!!!

## Compatible Applications

Flow Cytometry  
Immunocytochemistry  
Immunohistochemistry  
Immunofluorescence  
ELISA  
Cell based assays  
Immunopurification

**Secondary antibodies** can bind to primary antibodies or antibody fragments and offer significant advantages over directly labeled primary antibodies in target detection. By the sheer ability of multiple units of secondary antibody binding to a single primary antibody, they are capable of amplifying the target signal multifold and thereby increasing the detection sensitivity. BioVision is very pleased to offer an array of secondary antibodies for this very purpose. These secondary antibodies have been isolated from a variety of host species like **goat, rabbit, mice and sheep**. They are conjugated to commonly available probes like **HRP, FITC, TRITC, Cy3, Texas Red, AMCA, AP, Biotin**, etc. They can be used with colorimetric or fluorometric detection platforms.

These antibodies are compatible with numerous **detection, purification and sorting applications** as shown in the adjacent table. All these secondary antibodies are available in **convenient sizes and competitive prices**. Choose from the following list of secondary antibodies BioVision offers, to enhance your target signals.

## List of secondary antibodies from BioVision:

Secondary Antibody	Cat. No.	Size
Anti-Goat IgG, HRP-Linked Antibody	6403	500 µl
Anti-Mouse IgG, HRP-Linked Antibody	6402	500 µl
Anti-Rabbit IgG, HRP-Linked Antibody	6401	500 µl
Goat Anti-Human IgG (H&L) FITC	6915	100 µl
Goat Anti-Human IgG (H&L) HRP	6916	100 µl
Goat Anti-Mouse IgG (H&L) AMCA	6902	250 µl
Goat Anti-Mouse IgG (H&L) AP	6917	100 µl
Goat Anti-Mouse IgG (H&L) Biotin	6918	100 µl
Goat Anti-Mouse IgG (H&L) Cy3	6900	250 µl
Goat Anti-Mouse IgG (H&L) FITC	6919	100 µl
Goat Anti-Mouse IgG (H&L) HRP	6920	100 µl
Goat Anti-Mouse IgG (H&L) Texas Red	6901	250 µl
Goat Anti-Mouse IgG (H&L) TRITC	6921	100 µl

155 S. Milipitas Blvd, Milipitas, CA 95035

T: 408-493-1800 F: 408-493-1801

Toll Free: 800-891-9699 (US Only)

**BioVision**  
BioVision Incorporated



Secondary Antibody	Cat. No.	Size
Goat Anti-Mouse IgG (H&L) Unconjugated	6922	100 µl
Goat Anti-Rabbit IgG (H&L) AMCA	6905	250 µl
Goat Anti-Rabbit IgG (H&L) AP	6923	100 µl
Goat Anti-Rabbit IgG (H&L) Biotin	6924	100 µl
Goat Anti-Rabbit IgG (H&L) Cy3	6903	250 µl
Goat Anti-Rabbit IgG (H&L) FITC	6925	100 µl
Goat Anti-Rabbit IgG (H&L) TRITC	6926	100 µl
Goat Anti-Rabbit IgG (H&L) Unconjugated	6927	100 µl
Goat Anti-Rat IgG (H&L) Biotin	6910	250 µl
Goat Anti-Rat IgG (H&L) Cy3	6906	250 µl
Goat Anti-Rat IgG (H&L) FITC	6907	250 µl
Goat Anti-Rat IgG (H&L) HRP	6908	250 µl
Goat Anti-Rat IgG (H&L) TRITC	6909	250 µl
Mouse Anti-Goat IgG (H&L) AP	6928	100 µl
Mouse Anti-Goat IgG (H&L) FITC	6929	100 µl
Mouse Anti-Goat IgG (H&L) TRITC	6930	100 µl
Rabbit Anti-Goat IgG (H&L) Cy3	6911	250 µl
Rabbit Anti-Sheep IgG (H&L) FITC	6913	100 µl
Rabbit Anti-Sheep IgG (H&L) HRP	6912	100 µl
Rabbit Anti-Sheep IgG (H&L) TRITC	6914	100 µl
UltraPolymer Donkey Anti-Goat IgG (H&L) HRP	A1276	5 ml, 15 ml
UltraPolymer Donkey Anti-Sheep IgG (H&L) HRP	A1279	5 ml, 15 ml
UltraPolymer Goat Anti-Human IgG (H&L) HRP	A1278	5 ml, 15 ml
UltraPolymer Goat Anti-Mouse IgG (H&L) HRP	A1275	5 ml, 15 ml
UltraPolymer Goat Anti-Rabbit IgG (H&L) HRP	A1274	5 ml, 15 ml
UltraPolymer Goat Anti-Rabbit/Mouse IgG (H&L) HRP Cocktail	A1280	5 ml, 15 ml
UltraPolymer Goat Anti-Rat IgG (H&L) HRP	A1277	5 ml, 15 ml

Please visit [www.biovision.com](http://www.biovision.com) for a comprehensive overview on our products!